

EN

turn **all**® system

Automated Lateral Turning System




- **Automatic side positioning system 0-30° for pressure ulcer prevention or treatment**
- **Comprehensive pressure relief prevention of pressure injuries / ulcers treatment of existing injuries / ulcers**
- **Can be integrated into any care bed with removable mattress**
- **Time saving, results in more care time available for the patient**
- **Comfortable and doesn't interrupt normal sleep**
- **Side Support System slight lifting of the mattress, on the reverse side of the primary lift side, provides support and reassurance as well as decreases the risk of shear.**

LEVABO®



The Turn All® System turns a standard bed into an automated lateral turning bed



The new Levabo Turn All® lateral turning system is developed for acute and community settings .

The system is to support clinicians and carers in the regular turning of primary immobile patients, who are at the risk of developing pressure ulcers and for those who have existing pressure ulcers.

Clinical Efficiency

Clinical practice guidelines recommend regular repositioning of at least two hours for patients with reduced mobility (NICE GUIDELINE 2015).The 2019 International pressure ulcer guidelines (EPUAP,NPIAP, PPPIA) recommend the use of the 30 lateral side lying position.

Easy to Use

- Easy to set up and use
- On /off control button"
- 3 x turning cycles to choose from: 30/30, 60/60, Or 90/90minutes.
- Time setting to select the length of time that the automated lateral turning is required.
- CPR function for rapid deflation.
- The system gives an audible alarm if power is lost.



Control Box





The Turn All® System provides more time for caregiving.

Fits on every bed

Fits any standard bed where it is placed under the existing mattress.
Can be used with Dynamic mattresses and profiling beds.

Consistent and Controlled Turning

The system automatically turns the patient to gently and quietly reducing the need for manual handling intervention from care givers.



Side Support System

- Turn All® turning system is equipped with a support system for the opposite side of the primary lift. The support system eliminates the feeling of sliding away from the primary lift side, and reduces the risk of shear.



Benefits of the Turn All® system

- Comprehensive pressure relief
- Prevention of pressure injuries / ulcers
- Treatment of pressure injuries/ulcers),
- Comfortable and doesn't interrupt normal sleep
- Quiet when used
- Suitable for bariatric Management, and patients can be carefully and easily turned thereby protecting the health and wellbeing of staff working with them
- Gentle positioning even for patients in pain
- Slight lifting of the mattress, on the reverse side of the primary lift side, provides support and reassurance as well as decreases the risk of shear
- Time saving results in more care time available for the patient
- Improved microclimate and blood circulation
- Robust, durable and reliable
- Can be integrated into any care bed with removable mattress
- Provides consistent "head to toe" support
- Manual handling requirements are reduced.

Product items

Model	Size in cm.	Item no.	Max weight in Kg	Application
Turn All® System - Multi patient TPU 3C	85 x 200	70128	300	Set of 1 Turn All® lifting mattress mounted in Turn All® cover and 1 pump
Turn All® Inner lifting mattress Multi patient	85 x 200	70113	300	The active lifting part of the Turn All® System.
Turn All® Cover TPU	85 x 200	70114	NA	Cover for the Turn All® Mattress.
Turn All® Pump 3C	27 x 13 x 10	76117	NA	Pump for the Turn All® Mattress.

The "Turn All" system is supplied with a 5 year manufacturer's warranty on the pump and 2 year manufacturer's warranty on the mattress.



Distributor

Unit 7, Septimus Buildings
Hawkfield Business Park
Whitchurch
Bristol
BS14 0BN
TEL : +44 0117 9640818
EMAIL: sales@innovationrehab.co.uk
WEB: www.levabo.com

www.levabo.com
www.facebook.com/LevaboAps