
Effective Pressure Area Care BARIATRIC



EHOB® products have been clinically evaluated as part of effective protocol for the prevention and treatment of pressure injuries.



EHOB® has provided effective solutions for prevention and treatment of pressure injuries for more than 30 years.

We offer an extensive range of highly effective pressure relieving products to fit your specific needs.

EHOB® Bariatric Heel Protectors



4420A Heel elevator TPU Custom XL

- Complete heel elevation
- Fully adjustable
- Suitable for multipatient use
- Silicone Buttons
- Optional Extender Straps



TRUVUEXW060 TruVue Heel Protector XL

- Complete heel elevation
- Fully adjustable
- Machine washable
- Special mix of fibres to ensure optimal pressure redistribution
- Optional Antirotation Wedge

EHOB® Bariatric Product Range:

| Order number | Description | Size (inflated) | Weight Limit |
|-------------------------------|--|--|----------------------|
| 4010A | EHOB Bariatric Overlay Mattress | 192 cm x 120 cm x 8 cm | 360 kg |
| 4250A 4250B | EHOB Bariatric Cushion EHOB Bariatric Cushion with pump | 71 cm x 56 cm x 5 cm 71 cm x 56 cm x 5 cm | 320 kg 320 kg |
| 4420A 4420B TRUVUEXW060 | EHOB Heel Elevator TPU Custom XL As above including LAD pump TruVue Heel Protector - Extra Large | Calf Size 46 - 64 cm Calf Size 46 - 64 cm Calf Size 46 - 64 cm | N.A. N.A. N.A. |
| 4420AS 504MI | Extender Straps for 4420 Electric Inflation Pump for 4010 | | |

Can be reused following local policy guidelines for cleaning and disinfection. All EHOB products in the Static Air range incorporate an antimicrobial agent that acts as an effective barrier against the following microorganisms: Methicillin-Resistant Staphylococcus Aureus (MRSA); Clostridium difficile; Escherichia coli; Klebsiella pneumonia; Pseudomonas aeruginosa; Salmonella choleraesuis; Candida albicans; Enterococcus faecium; Aspergillus niger and Bacillus subtilis.



